



OrthoID Guidebook

Quick-start Guide to the Using the iOS App
and Participating in the Community

December 5, 2016



Introduction to OrthoID

Welcome to OrthoID

Doctors and industry experts use **OrthoID** to improve patient outcomes by quickly, accurately, and collaboratively identifying implanted orthopedic devices.


Contents

- | | |
|----------------------------|--------------------------|
| 1: Create an Account | 5: Make a Recommendation |
| 2: Download the App | 6: Join the Discussion |
| 3: Browse the Case Library | 7: Monitor Your Activity |
| 4: Create a Case | 8: Manage Settings |

1: Create an OrthoID Account

In order to begin using **OrthoID**, you'll first need to create an account:

- Navigate to <https://www.ortho.id> and select "Join the Community"
- Complete the registration form, making sure to enter the invitation code
- Select the "Join" button to confirm



NOTE: Be sure to enter your Invitation Code correctly. If you do not have an Invitation Code, you can create an account but will be unable to login to the OrthoID app until your account has been activated by the OrthoID team.

2: Download the App


After creating an account, a confirmation email will be sent to the email address that was entered on the registration form:

- Select the link in your confirmation email
- Download **OrthoID** from the AppStore
- Login to **OrthoID** with the email and password you registered



3: Browse the Case Library

On the 'Cases' navigation screen, the **OrthoID** case library is fully browseable. Each case is represented by a thumbnail view of the first case photo, a case number, the status of the case (open/-closed) and community's recommendation about the implanted device.




The **OrthoID** Case Library can be filtered by joint type, device manufacturer, date of case creation and more. Tap the orange 'Filter' button to open the filter selection screen and select a filter type. To search for a specific term, simply tap the Search bar and input your search terms

4: Create a Case

The first step to create a case is to select the 'Create' navigation screen, then:

- Select the a joint (Knee, Hip, Elbow, Ankle or Shoulder)
- Take a photo of the X-Ray, or select an X-Ray photo from your photo album
- Agree to the photo upload guidelines
- Complete & save the case




NOTE: Photos found to contain objectionable or non-medical content, disclose PHI (personal health information) or otherwise misuse the OrthoID system may be deleted without warning.

5: Make a Recommendation

Near the top of the 'Recommendation' tab on a selected case, you can view case photos as well as the case profile. All of the community's recommendations are available for voting toward the lower section of the 'Recommendation' tab, in addition to a section that allows you to submit your own recommendation.




Scroll left and right through photo thumbnails or tap to enlarge



Tap a checkmark to vote for an existing recommendation.

Tap the + and enter a Manufacturer, Model, and your comments, then tap Save to Submit a New ID.



SUBMIT A NEW ID +

MANUFACTURER
Blue

Yellow
Orange
Green
Blue
Purple
Pink
Black

MODEL
Unknown


EXPLAIN YOUR RECOMMENDATION
Tell us why you recommended this ID

Save

6: Join the Discussion

Community comments and contributions to a case are logged in the 'Discussion' tab on a selected case. Every vote, recommendation and comment on a case is marked with the contributor's name, organization and a time/date stamp. This helps you to gauge reliability of community input and can help the case creator come to a final decision.


To add a comment to the discussion, tap the text field at the bottom of the screen, enter your message, and then tap the orange arrow to submit.



7: Monitor Your Activity

To help you track what has changed since you last opened the **OrthoID** app, the 'Activity' navigation screen presents a brief summary of recent activity on all cases you have created or participated in.


You can quickly check case status, leading community recommendations and how many new votes, recommendations or comments were made on a case. Tap a case preview to see more details.



8: Manage Settings

In the 'Settings' navigation screen, you can review a summary of your **OrthoID** contributions and manage your preferences for the Tip Toolbar and Notifications.

There are three notification settings: you can choose to receive notifications for activity on all cases, only cases you created, or no notifications at all.



Swipe up or down to scroll the menu



Tap to edit Tip Toolbar and Notification settings